

# TGX MEDICAL

## STERILIZER LOAD LIST - ACTIVE LOAD

LOT  
CONTROL



**Load Number:** 20060403-0009-02  
**Sterilizer:** Riverside-Steam 1  
**Load Creation Time:** 4/3/2006 9:30:59AM  
**Sterilization Start Time:** 4/3/2006 9:32:21AM

**Load ID:** 279  
**Sterilization Process:** STEAM 270 - 5 Minutes  
**User ID:** agold

CONTENTS				
Qty	Item Description	SET ID	BuildNumber	Adding User ID
	HIP RETRACTORS HOLMAN	R5555150-2	22549	agold
	Holder, Lens Biom	R5502400-1	22701	agold
	KNIFE-LIBERATOR	R5555700-1	22715	agold
	KNIFE-LIBERATOR	R5555700-2	22564	agold
	HEMOCLIP SET	R5252875-2	22901	agold
	LAP CHOLE EXTRAS	R1253075-8	22596	agold
	Lens, Biom Large	R5502900-3	22714	agold
	Lens, Nikon 20D	R5503600-2	22703	agold
	ORTHO MAJOR INSTRUMENTS	R1556850-2	22632	agold
	ORTHO MINOR INSTRUMENTS	R1556900-3	22547	agold
	Lens, Biom One-Hole	R5503100-1	22699	agold
	Lens, Nikon 20D	R5503600-3	22710	agold

**12:ITEMS IN LOAD**

**STEAM LOAD PARAMETERS**

Prevac: Chamber Pressure = 28-32P\*, Temp=270° Minimum, Sterilize = 5 min.  
 Gravity: Chamber Pressure = 17-22P\*, Temp=250° Minimum, Sterilize = 40 min.  
 Bowie Dick: Chamber Pressure = 28-32P\*, Temp=270° Minimum, Sterilize = 3.5 min.  
 \*(1st "S" after the "C's")  
 Place your initial next to each parameter on the printout.

**STERRAD LOAD PARAMETERS**

Did the <b>chemical indicator</b> and/or sterilization tape on all items turn <b>yellow</b> ?	YES	NO
Is there a <b>Process Complete</b> printed on the tape?	YES	NO
Are <b>all letters</b> on the printer tape colored <b>black</b> ?	YES	NO

If No, please comment \_\_\_\_\_

**STERIS LOAD PARAMETERS**

<b>DIAGNOSTIC CYCLE:</b>	PASS	FAIL
Routine Cycle: Temperature _____		( 50°- 56° )
Concentration _____		(175 or greater)
Exposure Time _____		(12 minutes)
Correct Chemical Indicator Color Change	YES	NO
Quick Connect Properly Attached at End of Cycle	YES	NO

**ETHYLENE OXIDE PARAMETERS**

Total Cycle Time: \_\_\_\_\_ (15.00 + hours)  
 Aeration (Time) : \_\_\_\_\_ (12.00 + hours)  
 Sterilize (Time) : \_\_\_\_\_ (1.00 + hours)  
 Sterilize (Temp) : \_\_\_\_\_ (Approx. 130°)

Load Number: 20060403-0009-02

Load ID:

279

Sterilizer: Riverside-Steam 1

Sterilization Process: STEAM 270 - 5 Minutes

Load Creation Time: 4/3/2006 9:30:59AM

User ID: agold

Sterilization Start Time: 4/3/2006 9:32:21AM

Bowie Dick is:   ACCEPTABLE           UNACCEPTABLE  
SIGNATURE \_\_\_\_\_

If UNACCEPTABLE, state why and the action taken \_\_\_\_\_

---

Load is:           ACCEPTABLE           UNACCEPTABLE

